



Ireland's Leading Sourcing Specialists...

OHAUS GUARDIAN 7000 HOTPLATE STIRRERS

Durable Hotplate Stirrers with Superior Safety Features and Intelligent Performance

OHAUS Guardian 7000 Hotplate Stirrers are designed with outstanding safety features and superior heating and mixing performance. The SmartHousing™ is chemical resistant with an innovative, easy-to-clean design that channels spills away from internal components and the angled glass control panel. The bright LCD display includes intuitive icons, temperature, speed and time settings - all designed to monitor performance from across the lab.



GUARDIAN 7000 HOTPLATE STIRRERS FEATURES

Protect your lab with patented SmartPresence™ & SmartLink™ technologies. These exclusive features automatically shut off the heater if no one is detected after a user set "time out" period has passed.

An industry-leading early detection system that uses two independent safety controls to continuously monitor the electronics and will shut off heating before an overtemperature condition occurs.

Precise stirring from 60 to 1600 rpm is made possible by the design of a powerful motor, strong magnet and software controlled ramp rate, ensuring secure magnetic coupling for viscous applications.

Call us today for more information or a quotation.