

Contact Details

Email: niall.fitzpatrick@fitzpatricklabsupplies.ie

Tel: +353 (0) 87 343 1487

Email: eamonn.carey@fitzpatricklabsupplies.ie

Tel: +353 (0) 87 738 5011

PPE - Personal Protective Equipment

PE Apron

Product details

- extra strong
- embossed surface
- neck loop
- PE-ribbon ties in waist

Material

- polyethylene
- moisture-resistant, solid
- $\approx 0.05 \text{ mm } (\pm 5 \%)$

Norms & guidelines

- CE acc. to MDD – Directive 93/42/ECC + 2007/47/EC – Class 1 product
- disposable product
- latex-free
- non sterile
- suitable for contact with food

Colour	Size	Order Code	Pack Size
White	80 x 130cm	20100703	500



Filtering Facepiece Cobra Foldy FFP3/V D

Product information:

- Construction of mask in 3D panel design
- Flexible nose clip
- With comfortable, continuously adjustable, cycled harness and unique Comfort-Harness clip
- Metal-free band holder-strap
- Specially designed for wearer of EKAStU spectacles and goggles
- Individually packed to avoid contamination during storage
- Protection against toxic and highly toxic particles up to 50 times the NPF value
- Fulfilling the requirements against high dust pollution and therefore are authorized wearing the marking D.
- Protection class DIN EN 149:2001 FFP3 NR D & (EU) 2016/425
- In original packaging: 5 years shelf life
- **Ordering Code 419 283 - 12 pieces per pack.**



PPE - Personal Protective Equipment

Re-Usable Face Mask

Outer layers:

3 layer mask
100% Cotton, Single Jersey, 180g/m2, Okotex 100 Class1

Middle layer:

100% Polyester, Polar fleece, 150g/m2, Okkotex 100 Class1

Attachment :

Polyester elastics role to attach around the ear

Made under SA800/BSCI to ensure a safe and secure environment

Manufactured under ISO 9001:2015 and in compliance with all hygienic guidelines

Made from sustainable materials through ecologically friendly production within ISO1401 environmental standard factories

Size	Colour	Pack Size
19.5cm x 16cm	Black	5 per pack



Face Shield

Face Shield Description:

The Face Shield is composed by the support, the full faced shield and the elastic headband.

Face Shield Support:

Material - Polypropylene
Weight - 22,5 ± 0,7 g
Colour - Black

Full Face Shield:

Material - PET
Weight - 32,1 ± 1,0 g





PPE - Personal Protective Equipment

Disposable Isolation Gown

(Non-Medical).

Description

45g PP Non-Woven with PE Film

Dimension

120 cm x 140 cm

Scope of Application

This product is used as an ordinary isolation gown for outpatient and inspection rooms of medical institutions.

Instructions

- Necklink with velcro, easy to close the clothes.
- One piece design for easy putting on and off.
- Tight cuffs,the elasticity is tightened.
- Waist design and the lacing method makes the clothes more fit to the body.

Packaging

One piece in a polybag and vacuum packing
100 pieces in a carton



Fitzpatrick Laboratory Supplies

Your Sourcing Specialists

Bespoke, flexible, efficient ... it's what we do.

Contact us today to discuss all your Laboratory needs.