



accu-jet® pro Pipette Controllers

All with one hand: select the delivery mode (gravity-delivery / blow-out) and adjust the motor speed range with the thumb; use variable button pressure for finest control of filling and delivery speed. The ergonomic handgrip weights only 190 g, perfectly balanced design - for fatigue-free pipetting even in prolonged operations.

Supplied with two battery compartment covers. A second cover to place the instrument in an inverse position on the bench, to avoid contamination of the pipette adapter.

- Ergonomic design.
- Adjustable motor speed limiter - a slower speed setting allows even more sensitive control of filling and dispensing rates for lower volume pipettes.
- A flashing LED indicates recharging of the storage battery.
- Approx. 8 h of continuous pipetting (with a 10 ml pipette) without recharging.
- Direct outlet of liquid vapours protects against corrosion to ensure long instrument life.
- A 0.2 µm hydrophobic membrane filter and a built-in check valve protect against liquid penetrating.
- Perfect fit for glass and plastic pipettes from 0.1 to 200 ml.

Description	Colour	Order Number
accu-jet® pro	Dark Blue	5281090
accu-jet® pro	Magenta	5281010
accu-jet® pro	Green	5281015
accu-jet® pro	Royal Blue	5281020
accu-jet® pro, without recharger	Dark Blue	5281025

Call Fitzpatrick Laboratory Supplies today for more information, our experts are on hand to discuss all your requirements.

Follow our LinkedIn page for new product and industry news.

 <https://www.linkedin.com/company/fitzpatrick-laboratory-supplies>

niall.fitzpatrick@fitzpatricklabsupplies.ie

Tel: +353 (0) 87 343 1487

eamonn.carey@fitzpatricklabsupplies.ie

Tel: +353 (0) 87 738 5011