

DWK Life Sciences WHEATON® COMPLETEPAK™



DESIGNED BY YOU.
CRAFTED BY US.



DURAN
WHEATON
KIMBLE

Excellence in your hands

WHEATON® COMPLETEPAK

The WHEATON® COMPLETEPAK product line is a range of crimp top vials, rubber stoppers and seals that are ready-to-use and are off-the-shelf, which can be customized for any application. This reduces our customer's supply chain due to our one source solution capability. Each kit comes with a specific United States Pharmacopeia (USP) & European Pharmacopoeia (EP) certificate showing that the product meets or exceeds critical USP / EP standards that are enforced by the Federal Drug Administration (FDA) & European Directorate for the Quality of Medicines and Healthcare (EDQM). This kind of solution allows for our customers to focus on their core competency which is manufacturing, compounding and developing drug products.

Target Applications

- Compounding pharmaceuticals
- Aseptic filling of injectable drug products
- Drug discovery
- Drug research
- Freeze drying / lyophilized drug products
- Large and small molecule drug packaging

Key Benefits

- Sterilized, off-the-shelf, ready-to-use (RtU) products
- Variety of Kitted and stand-alone sterile components — a solution for any application
- Kitted solutions feature equal amounts of RtU components for reduced scrap rates
- Complete manufacturing & certificate of analysis traceability
- USP <788>/EP 2.9.19 Particulate Matter in Injections, USP <85>/EP 2.6.1, The Bacterial Endotoxins Test and USP <71>/EP 2.6.14 The Sterility Test, USP <660>/EP 3.2.1 Glass Containers for Pharmaceutical Use, and USP <381>/EP 3.2.9 Elastomeric Closures used in Parenteral Preparations
- CompletePAK Amber for light sensitive content packaging
- CompletePAK ModPAK Bulk configurations
- Ready-to-Use, Just-in-Time, minimum 3 month shelf life from time of shipment
- Container Closure Integrity (CCI) certification available upon request

Materials

- Vials are designed for parenteral applications and manufactured from low extractable borosilicate glass
- All stoppers are manufactured from an ultrapure bromobutyl latex free formulation with extremely high chemical purity and low gas permeability
- Omniflex stoppers offer low levels of residual binding to maximize extraction of active pharmaceutical ingredients
- Omniflex stopper elasticity delivers a more effective barrier to optimize container closure integrity and content protection
- Igloo lyophilized stoppers are suited for biologics and content with sensitivities to silicon



ENHANCING THE EFFICIENCY OF DRUG DISCOVERY AND DRUG COMPOUNDING



TOP REASONS WHY YOU SHOULD CHOOSE COMPLETEPAK FOR YOUR APPLICATION:

1. Products supplied ready-to-use to improve workflow
2. Single source options reduce lead times & Simplify administrative effort
3. Focuses your efforts on critical core competencies
4. Eliminates capital investments in sterile processing equipment
5. Frees up valuable clean room space
6. Validated processes & Certificate of Analysis Traceability assures quality production outcomes
7. Saves valuable time and money while maintaining a high quality product and streamlining process

COMPLETEPAK Kits

Clear Kit Cat. No.	Amber Kit Cat. No.	Description (pieces/kit)	Qty / Case	ModPAK Cat. No. (Clear Only)	Cases / ModPAK
W990000-CP	W990000-CPA	2mL Sterile Vial (210), 13mm Ultra Pure Sterile Serum Stopper (230), Sterile Blue Seal (230)	1	W990000-CPMP	6
W990001-CP	W990001-CPA	2mL Sterile Vial (210), 13mm Omni Flex 3G Sterile Serum Stopper (230), Sterile Blue Seal (230)	1	W990001-CPMP	6
W990002-CP	W990002-CPA	2mL Sterile Vial (210), 13mm OmniFlex 3G Sterile Igloo Lyo Stopper (230), Sterile Blue Seal (230)	1	W990002-CPMP	6
W990003-CP	W990003-CPA	2mL Sterile Vial (210), 13mm Ultra Pure Sterile Serum Stopper (230), Sterile Red Seal (230)	1	W990003-CPMP	6
W990004-CP	W990004-CPA	2mL Sterile Vial (210), 13mm Omni Flex 3G Sterile Serum Stopper (230), Sterile Red Seal (230)	1	W990004-CPMP	6
W990005-CP	W990005-CPA	2mL Sterile Vial (210), 13mm OmniFlex 3G Sterile Igloo Lyo Stopper (230), Sterile Red Seal (230)	1	W990005-CPMP	6
W990006-CP	W990006-CPA	5mL Sterile Vial (200), 20mm Ultra Pure Sterile Serum Stopper (220), Sterile Blue Seal (220)	1	W990006-CPMP	6
W990007-CP	W990007-CPA	5mL Sterile Vial (200), 20mm OmniFlex 3G Sterile Serum Stopper (220), Sterile Blue Seal (220)	1	W990007-CPMP	6
W990008-CP	W990008-CPA	5mL Sterile Vial (200), 20mm OmniFlex 3G Igloo Lyo Stopper (220), Sterile Blue Seal (220)	1	W990008-CPMP	6
W990009-CP	W990009-CPA	5mL Sterile Vial (200), 20mm Ultra Pure Sterile Serum Stopper (220), Sterile Red Seal (220)	1	W990009-CPMP	6
W990010-CP	W990010-CPA	5mL Sterile Vial (200), 20mm OmniFlex 3G Sterile Serum Stopper (220), Sterile Red Seal (220)	1	W990010-CPMP	6
W990011-CP	W990011-CPA	5mL Sterile Vial (200), 20mm OmniFlex 3G Sterile Igloo Lyo Stopper (220), Sterile Red Seal (220)	1	W990011-CPMP	6
W990012-CP	W990012-CPA	10mL Sterile Vial (158), 20mm Ultra Pure Sterile Serum Stopper (220), Sterile Blue Seal (220)	1	W990012-CPMP	6
W990013-CP	W990013-CPA	10mL Sterile Vial (158), 20mm OmniFlex 3G Sterile Serum Stopper (220), Sterile Blue Seal (220)	1	W990013-CPMP	6
W990014-CP	W990014-CPA	10mL Sterile Vial (158), 20mm OmniFlex 3G Sterile Igloo Lyo Stopper (220), Sterile Blue Seal (220)	1	W990014-CPMP	6
W990015-CP	W990015-CPA	10mL Sterile Vial (158), 20mm Ultra Pure Sterile Serum Stopper (220), Sterile Red Seal (220)	1	W990015-CPMP	6
W990016-CP	W990016-CPA	10mL Sterile Vial (158), 20mm OmniFlex 3G Sterile Serum Stopper (220), Sterile Red Seal (220)	1	W990016-CPMP	6
W990017-CP	W990017-CPA	10mL Sterile Vial (158), 20mm OmniFlex 3G Sterile Igloo Lyo Stopper (220), Sterile Red Seal (220)	1	W990017-CPMP	6

COMPLETEPAK Sterile Components*

Cat. No.	Description	Qty / Case
Clear Vials		
W986000-S	2mL Clear Sterile Serum Vial	210
W986001-S	5mL Clear Sterile Serum Vial	200
W986002-S	10mL Clear Sterile Serum Vial	158
Amber Vials		
W986000A-S	2mL Amber Sterile Serum Vial	210
W986001A-S	5mL Amber Sterile Serum Vial	200
W986002A-S	10mL Amber Sterile Serum Vial	158

Cat. No.	Description	Qty / Case
Stoppers		
W986003-S	13mm OmniFlex 3G Sterile Serum Stopper	230
W986004-S	13mm Ultra Pure Siliconized Sterile Serum Stopper	230
W986005-S	13mm OmniFlex 3G Sterile Igloo Lyo Stopper	230
W986006-S	20mm OmniFlex 3G Sterile Serum Stopper	220
W986007-S	20mm Ultra Pure Siliconized Sterile Serum Stopper	220
W986008-S	20mm OmniFlex 3G Sterile Igloo Lyo Stopper	220
Seals		
W986009-S	13mm Sterile Flip Red Seal	230
W986010-S	13mm Sterile Flip Blue Seal	230
W986011-S	20mm Sterile Flip Red Seal	220
W986012-S	20mm Sterile Flip Blue Seal	220

*Sterile components are available in quantities of 2 or less. For quantities of 3 or more, contact us directly.



WHEATON® is a registered trademark of DWK Life Sciences LLC. COMPLETEPAK is a trademark of DWK Life Sciences LLC.
Copyright © 2019 DWK Life Sciences LLC. WBRO_058 06/2020

DWK Life Sciences GmbH
Hattenbergstrasse 10
55122 Mainz
Germany

Phone: +49 (0)6131/1445 4131
Fax: +49 (0)6131/1445 4016
sales@dwk.com
www.dwk.com