

Technical Specification

PRODUCT : Sodium acetate trihydrate RE - Pure**CODE : 366100****MOLECULAR WEIGHT :136,08****CAS N° : 6131-90-4****FORMULA : CH₃COONa.3H₂O**

TEST	U.M.	SPECIFICATION
Description	-	White or almost white,cryst. powder or colourless cryst.
Identification (I.R.)	-	Positive
Appearance of solution	-	Conform
Assay (on dry)	%	99.0 - 101.0

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in