



Sodium acetate trihydrate

Synonyms:

Acetic acid sodium salt trihydrate

CH₃COONa.3H₂O

Classification

Danger

Molecular weight **136.08** ONU

CAS **6131-90-4** Transport hazard class:

EEC-N: **612-115-9** Packing

Sodium acetate trihydrate

> ERBAPharm - According to pharmacopoeia :
Ph.Eur.-USP-FU-BP

Pharmaceutical

Description	White or almost white,cryst. powder	-
.	or colourless cryst.	-
Identification (I.R.)	Positive	-
Appearance of solution	Conform*	Ph.Eur.
Reducing substances	Conform*	Ph.Eur.
K	Conform*	USP-NF



Sodium acetate trihydrate> ERBAPharm - According to pharmacopoeia :
Ph.Eur.-USP-FU-BP

Pharmaceutical

pH sol. 5% at 25°C	7.5 - 9.0	-
Loss on drying (130°C)	39.0 - 40.5	%
Chloride	<=200*	ppm
Heavy metals (as Pb)	<=10*	ppm
Sulfate	<=50*	ppm
Not soluble matter	<=500*	ppm
As	<=2*	ppm
Calcium + Magnesium	<=50*	ppm
Fe	<=10*	ppm
Al	<=0.2*	ppm
Assay (on dry)	99.0 - 101.0	%

ID	Size	Packaging
000000000000366207	1 kg	Plastic bottle
000000000000366209	5 kg	Plastic tank
000000000000366205	25 kg	Plastic bucket
000000000000366206	25 kg	Plastic bucket
000000000000366204	50 kg	Fibre drum



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA Reagents
operates with a Certified
Quality Management
System